

EXCEED[®] 500_{xr}

CUSTOM LIQUID

EXCEED[®] 500_{xr} CUSTOM LIQUID FOR SOYBEAN

This inoculant is intended to be applied on-seed. It is a culture of highly effective strain of rhizobia (nitrogen-fixing bacteria) that will colonize the plant roots to gather and fix the “free” nitrogen of the air, making it available to the plants.

Product 1015

Treats: 25,000# (11,340 kg) of seed or 500 units

Contents: One 12.0 liter (405.77 fl oz) Bladder bag of inoculant and One 2.8 liter (94.6 fl oz) Bladder bag of extender

Package Configuration: 1 package per case, 49 cases per pallet

Product 1016

Treats: 2,500# (1,134 kg) of seed or 50 units

Contents: Four 1.15 liter (38.88 fl oz) Bladder bag of inoculant and four 0.35 liter (11.83 fl oz) Bladder bags of extender

Package Configuration: 4 sets of one inoculant bladder and one extender bladder per case, 54 cases per pallet

**Application rates are based upon 2,800 seeds per pound (0.45 kg). **

- Prepare only as much Exceed[®] 500_{xr} Custom Liquid for Soybean as will be used that day.
- Consult with a Visjon Biologics representative for a complete list of compatible seed treatments.

The product density for Exceed[®] 500_{xr} Custom Liquid for Soybean is 8.6 lbs/gal (1.03 g/mL)

Storage of Exceed[®] 500_{xr} Custom Liquid for Soybean

- Always read and follow instructions.
- Store in a cool, dry place out of direct sunlight.
- Store product between 40° and 77° F (4° and 25° C.)
- Open package only when ready to use.
- Use before expiration date.
- Do not store in an area which is to be fumigated.
- Small amounts of product which cannot be used according to instructions may be disposed of in a landfill.

Expiration Date:

Exceed[®] 500_{xr} Custom Liquid for Soybean carries a one-season shelf life and expires the date stamped on the case.

Information regarding the contents and levels of metals in this product is available on the internet at:
<http://www.aapfco.org/metals.html>

Limited Warranty:

VISJON BIOLOGICS GUARANTEES THIS CULTURE TO PRODUCE SATISFACTORY NODULE FORMATION AND NITROGEN FIXATION UNDER FAVORABLE CONDITIONS OF SOIL AND CLIMATE WHEN USED BEFORE EXPIRATION DATE AND APPLIED UNDER THE MANUFACTURER'S SPECIFICATIONS. OR PURCHASE PRICE WILL BE REFUNDED. A VISJON BIOLOGICS REPRESENTATIVE MUST BE NOTIFIED OF ANY FIELD COMPLAINT WITHIN FIFTY (50) DAYS AFTER PLANTING. MANY FACTORS OTHER THAN NODULATION AFFECT CROP PERFORMANCE. VISJON BIOLOGICS ASSUMES NO RESPONSIBILITY FOR LOSS OR PARTIAL LOSS OF CROP FROM ANY CAUSE WHATSOEVER. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. THE LIMITED WARRANTY IS VOID WHERE PROHIBITED BY LAW. IF YOU HAVE ANY QUESTIONS ABOUT THIS PRODUCT OR ANY OTHER VISJON BIOLOGICS PRODUCT, CALL 940-264-0343 FOR ASSISTANCE.

CONTAINS NONPLANT FOOD INGREDIENTS

Minimum Guaranteed Analysis:

Active Ingredient: <1% by Weight

Bradyrhizobium japonicum per ml: 1x10¹⁰ cfu/mL

Inactive Ingredient: >99% by Weight

KEEP OUT OF REACH OF CHILDREN

This product does not meet the regulatory definition of a hazardous chemical under 29 CFR 1910.1200. However, good industrial hygiene practices should be used in the manufacture and handling of the product.

Precautionary Information:

While Exceed[®] 500_{xr} Custom Liquid for Soybean is non-toxic it is recommended that the applicator avoid inhalation. Visjon Biologics recommends the use of gloves and eye protection while using Exceed[®] 500_{xr} Custom Liquid for Soybean.

This product is designed to be custom applied through an approved application device by a retail location that has been trained in the proper methods of handling and applying Exceed[®] 500_{xr} Custom Liquid to soybean seed. Please read all directions and specifications carefully before using this product. This product contains live organisms that are sensitive to chemicals and harsh environmental conditions.

FOR RETAIL APPLICATION ONLY.

Directions for Application:

- Ensure that the treater has been properly calibrated.
- Shake each container well before use.
- Pour each container into mix tank.
- Agitate or recirculate liquid mixture to ensure proper distribution of bacteria and extender.
- Apply at a rate of 2.0 fl oz (60 mL) per cwt or 1.0 fl oz (30 mL) per 140,000 seeds.



1092 S Graham ST
Henrietta TX 76365
940-264-0343 | 940-264-0344 fax
www.visjonbiologics.com